

MAA519Rb21

Monoclonal Antibody to Apolipoprotein A1 (APOA1)

Organism Species: *Oryctolagus cuniculus* (Rabbit)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]

Immunogen: APOA1, Rabbit

Clonality: Monoclonal

Clone number: H6

Host: Mouse


Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 500µg/mL

UOM: 200µg



Sample: Recombinant APOA1, Rabbit

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant APOA1 (Asp25~Gln265) expressed in *E.coli*.

Accession No.: RPA519Rb01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHSGS-DEPRSS WDKIKDFATV YVDTVKDSGR EYVAQFEASA FGKQLNLKLL
DNWDSLSTV SKLQEQLGPV TQEFWDNLEK ETEGLREEMN KDLQEVQRKV
QPFLDEFQKK WQEEVERYRQ KVEPLGAELG ESARQKLTEL QEKLSPLAEE
LRDSARTHVD TLRTKLAPYS NELQRLAARL ESIKEGGGAK LAEYQAKARE HLSVLSEKAR
PALEDLRQGL LPVLESFKAS VQNVVDEATK KLNTQ

[ANTIBODY SPECIFICITY]

The antibody is a mouse monoclonal antibody raised against APOA1. It has been selected for its ability to recognize APOA1 in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.