




A93353Hu01
Polyclonal Antibody to Slit Homolog 3 (Slit3)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

4th Edition (Revised in September, 2012)

[PRODUCT INFORMATION]

Immunogen: Slit3
Clonality: Polyclonal
Host: Rabbit
Species Reactivity: Human
Ig type: Rabbit IgG
Purification: Antibodies are purified by target protein affinity chromatography.
Applications: WB, ICC, IHC-P, IHC-F
Form: Liquid
Size: 100µg



Sample: Recombinant human Slit3

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant human Slit3 (Lys1405~Glu1512) expressed in *E.coli*.
Molecular Weight: 17.6 kDa
USCN accession No.: P93353Hu01
Sequence: The target protein is fused with two N-terminal Tags, His-tag and S-tag and its sequence is listed below.
MHHHHHHSSGLVPRGSGMKETAAAKFERQHMDSPDLGTDDDDKAMADIGSEF-KNDSAN ACSAFKCHHG QCHISDQGEP
YCLCQPGFSG EHCQQENPCL GQVVREVIRR QKGYASCATA SKVPIMECRG GCGPQQCQPT RSKRRKYVFQ CTDGSSFVEE VE



[ANTIBODY SPECIFICITY]

Anti Slit3 is a rabbit polyclonal antibody raised against human Slit3. This antibody has been selected for its ability to recognize human Slit3 in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH 7.4, containing 0.02%Na₃N, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.

