

A94771Mu01
Polyclonal Antibody to Interleukin 11 Receptor Alpha (IL11Ra)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[PRODUCT INFORMATION]

Immunogen: IL11Ra

Clonality: Polyclonal

Host: Rabbit

Species Reactivity: Mouse


Ig type: Rabbit IgG

Purification: Antibodies are purified by target protein affinity chromatography.

Applications: WB, ICC, IHC-P, IHC-F

Form: Liquid

Size: 100µg



Sample: Recombinant mouse IL11Ra

[IMMUNOGEN INFORMATION]




Immunogen: Recombinant mouse IL11Ra (Trp59~Arg217) expressed in *E.coli*.

Molecular Weight: 18.9kDa

USCN accession No.: P94771Mu01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHSGSEF- WF RDGDSRLQG PDSGLGHLV LAQVDSPEG TYVCQTLDGV
SGGMVTLKLG FPPARPEVSC QAVDYENFSC TWSPGQVSGP PTRYLTSYRK KTLPGAESQR
ESPSTGPWPC PQDPLEASRC VVHGAEFWSE YRINVTEVNP LGASTCLLDV RLQSILR

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

[ANTIBODY SPECIFICITY]

Anti IL11Ra is a rabbit polyclonal antibody raised against mouse IL11Ra. This antibody has been selected for its ability to recognize mouse IL11Ra in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH 7.4, containing 0.02%Na₃N, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant mouse IL11Ra (Trp59~Arg217) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading buffer: 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na₃N 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.