

PAK207Bo01 **Polyclonal Antibody**

To A Disintegrin And Metalloproteinase With Thrombospondin 2 (ADAMTS2) Organism Species: Bos taurus; Bovine (Cattle)

Instruction manual

←17kDa

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

Western Blot

kDa 75

50

37

25

20

15

10

[PRODUCT INFORMATION]

Immunogen: ADAMTS2, Bovine

Clonality: Polyclonal

Host: Rabbit

Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 200µg/mL

Sample: Recombinant ADAMTS2, Bovine

UOM: 100µg [IMMUNOGEN INFORMATION]

Immunogen: Recombinant ADAMTS2 (Gln474~Pro610) expressed in E.coli.

USCN Accession No.: RPK207Bo01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence

is listed below.

MGHHHHHHSGS-QLPGLHY SMNEQCRFDF GLGYMMCTAF RTFDPCKQLW CSHPDNPYFC KTKKGPPLDG TMCAPGKHCF KGHCIWLTPD ILKRDGNWGA WSPFGSCSRT CGTGVKFRTR QCDNPHPANG GRTCSGLAYD FQLCNSQDCP



[ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against ADAMTS2. It has been selected for its ability to recognize ADAMTS2 in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200 Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant ADAMTS2 (Gln474~Pro610) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemilumescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN₃ 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.