

**PAP657Hu01**

**Polyclonal Antibody to Suppression Of Tumorigenicity 14 (ST14)**

**Organism Species: Homo sapiens (Human)**

***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

## **[ PRODUCT INFORMATION ]**

**Immunogen:** ST14, Human

**Clonality:** Polyclonal

**Host:** Rabbit


**Immunoglobulin Type:** IgG

**Purification:** Affinity Chromatography.

**Applications:** WB, ICC, IHC-P, IHC-F, ELISA

**Concentration:** 200µg/mL

**UOM:** 50µg



*Sample: Recombinant ST14, Human*

## **[ IMMUNOGEN INFORMATION ]**

**Immunogen:** Recombinant ST14 (Arg607~Asn852) expressed in *E.coli*.

**Accession No.:** RPP657Hu01

**Sequence:** The target protein is fused with two N-terminal Tags, His-tag and S-tag and its sequence is listed below.

MHHHHHHSSG LVPRGSGMKE TAAAKFERQH MDSPDLGTDD DDKAMADIGS EF-  
RSFTRQARVV GGTDADDEGEW PWQVSLHALG QGHICGASLI SPNWLVSAAH  
CYIDDRGFY SDPTQWTAFL GLHDQSQRSA PGVQERRLKR IISHPPFNDF TFDYDIALLE  
LEKPAEYSSM VRPICLPDAS HVFPAGKAIW VTGWGHTQYG GTGALILQKG EIRVINQTTT  
ENLLPQQITP RMMC VGFLSG GVDSCQGD SG GPLSSVEADG RIFQAGV VSW  
GDGCAQRNKP G VYTRLPLFR DWIKEN

## **[ ANTIBODY SPECIFICITY ]**

The antibody is a rabbit polyclonal antibody raised against ST14. It has been selected for its ability to recognize ST14 in immunohistochemical staining and western blotting.

## **[ APPLICATIONS ]**

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

## **[ CONTENTS ]**

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02% NaN<sub>3</sub>, 50% glycerol.

## **[ STORAGE ]**

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.