


300927\_H



Current Data Parameters  
NAME 20150211  
EXPNO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 20150211  
Time 14.24  
INSTRUM spect  
PROBHD 5 mm TXI 1H/D-  
PULPROG zg30  
TD 65536  
SOLVENT CDCl3  
NS 16  
DS 2  
SWH 10330.578 Hz  
FIDRES 0.157632 Hz  
AQ 3.1720407 sec  
RG 575  
DW 48.400 usec  
DE 6.00 usec  
TE 298.2 K  
D1 1.00000000 sec  
TD0 1

===== CHANNEL f1 =====  
NUC1 1H  
P1 8.85 usec  
PL1 1.00 dB  
SFO1 500.2530892 MHz

F2 - Processing parameters  
SI 32768  
SF 500.2500153 MHz  
WDW EM  
SSB 0  
LB 0.30 Hz  
GB 0  
PC 8.00

